

TV-DEX

CONVENIENT STEREO SOUND



- 1. TV-Controller
- 2. TV-Base
- 3. Transmission on/off
- 4. Volume up
- 5. Volume down
- 6. Room Off

WWW.WIDEX.COM



Printed by KLS/2016-05
9 502 1839 001 #01



WIDEX[®]
HIGH DEFINITION HEARING

TV-DEX

FOR ENJOYING TV



DEX[™]
A WIRELESS WORLD

A WHOLE NEW WORLD OF **WIRELESS** **SOUND**

The new DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.

They make communication between your hearing aids and your remote control, television or mobile phone easy.

TV-DEX

The TV-DEX is a user-friendly wireless assistive listening device specifically designed for enjoying TV and audio. You can enjoy up to 10 hours of uninterrupted, top quality TV or audio sound between recharges.

The main advantage of the TV-DEX is its real-time high-quality stereo sound. You experience TV or audio sound exactly as it was intended - without annoying distortion or more importantly - echoes.

Another unique feature of the TV-DEX is the 'Room Off' function. This allows you to temporarily switch the hearing aid's microphone off and hear the TV sound only - meaning you can enjoy your favourite TV programme without unnecessary background noise. Naturally, you can just as easily switch the 'Room Off' function off again if you prefer.

