

WIDEX CLEAR™ 440

STAY IN SYNC WITH LIFE



WIDEX[®]
HIGH DEFINITION HEARING

WIDEX CLEAR440-PA

A family of four is sitting on large, dark rocks by the ocean. On the left, a woman in a white long-sleeved shirt and white pants is laughing. In the center, a woman wearing sunglasses and a blue scarf is also laughing. To her right, a young girl with curly hair is smiling. On the far right, a man in a white short-sleeved shirt and khaki pants is sitting and talking to the others. The background shows the ocean and a rocky coastline under a bright sky.

NATURE
HAS GIVEN US
TWO EARS
DESIGNED
TO **WORK**
TOGETHER



A WORLD OF **NATURAL** SOUND

Nature has given us two ears designed to work together.

And like two ears, the new WIDEX CLEAR™440-PA hearing aids communicate with each other. Instantly and constantly.

The CLEAR440-PA model does this completely wirelessly. All functions in each hearing aid can be adjusted seamlessly. So you can experience clear sound as close to natural as possible.

FLEXIBLE

Being wireless, CLEAR440-PA makes it easy for you - via a uniquely designed range of accessories - to handle some of the more difficult listening situations such as watching television, listening to music or using your mobile phone.

CLEAR440-PA can be configured exactly the way you would like. So you can always hear as naturally as possible in different situations. Just as nature intended.

HEARING AID **MODEL** AND **COLOURS**

WIDEX CLEAR440-PA is available in a full range of up to fourteen fantastic colours!



Standard colours

	Warm beige		Midnight black
	Tan silk		Titan grey
	Cappuccino brown		Winter silver

Additional colours

	Summer gold		Shocking pink
	Copper brown		Lime green
	Pearl white		Sporty red
	Metallic blue		Mediterranean turquoise

NATURAL SOUND

At Widex we are only too aware that you want to hear as naturally as possible. Indeed, the high quality sound of our hearing aids is something we are both proud of and are known for.

With CLEAR440-PA, we have actively tried to take this a step further. By improving the way hearing aids communicate with each other, we have made sure that CLEAR440-PA sounds and works like no other hearing aid – as close to natural sound as possible.





NATURALLY **WIRELESS**

Having access to an array of advanced features in your hearing aid is one thing; to get the most out of them though, CLEAR440-PA uses new wireless technology, WidexLink. This guarantees that your hearing aids can communicate with each other instantly. **Just like two ears do.**

STAY IN TOUCH

Every time CLEAR440-PA hearing aids adjust to your listening environment, or every time you switch to a different program or different feature, they react without delay and **stay in contact.**

And if you lose your hearing aid or if the battery is exhausted, the Partner Monitor in CLEAR440-PA **alerts you immediately.**

Put simply, our wireless technology means real sound in real time.

EASY LISTENING

Even at the best of times, listening clearly can be a challenge when using your hearing aid. In particular, when going from one sound environment (such as soft sounds) to another (such as loud sounds).

CLEAR440-PA solves this challenge by making the change between listening environments – from soft sounds, to normal speech, to loud sounds – as smooth and comfortable as possible. This means **sounds are more detailed, natural and clear.**

CLEAR SOUND

Determining what is noise and what is speech has always been one of the biggest challenges facing hearing aid users. The IE Speech Enhancer in CLEAR440-PA lets you better understand speech in noisy situations and **gives speech the clarity it deserves.**

The visible part of the ear, the pinna, helps collect sound – a bit like a funnel. This means that the pinna helps you in determining which direction the sound is coming from. This can be difficult for hearing aid users. CLEAR440-PA features a unique Digital Pinna that **helps you locate sounds from any direction.**



FREEDOM OF CHOICE

A range of unique features give you the possibility of matching your hearing aids to your particular needs.

- FreeFocus is our answer to the challenge of hearing from different directions. With this feature you can **focus on sound in selected directions** – front, back, left and right – without having to turn your head.
- Our IE Zen **dedicated relaxation** and tone program plays random and harmonic tones in stereo. Whenever you feel like relaxing, you can hear them in a range of different tones or 'styles'.
- Your hearing care professional can help you adjust the tempo, volume and pitch of the tones to suit your taste.
- And with the SmartSpeak™ system, you are **informed verbally** with actual recorded speech, or tones if you prefer, of various hearing aid functions, such as low battery power and program mode.

So hearing naturally and clearly could not be easier.





THE COMPLETE HEARING SOLUTION

You can widen your listening experience by combining your CLEAR440-PA hearing aids with one of our assistive listening devices.

CLEAR WORLD OF SOUNDS

Watching television can sometimes be a challenge for hearing aid users, with unnecessary delay or echo. With our TV-DEX, high quality sound is transmitted with no delay straight to your CLEAR440-PA. Best of all, our wireless technology ensures that there is no echo so you can enjoy **uninterrupted sound** for up to ten hours between recharging the battery.

And you can listen with the hearing aid microphones switched off so surrounding sounds do not interfere with your listening experience.

The TV-DEX can also be used to hear music from your stereo or other audio source.

With the M-DEX, you can use your hearing aids as wireless headsets. They **connect directly** with any mobile phone, allowing you to talk and hear mobile phone conversations directly in your hearing aids.

And a unique feature of the M-DEX is what we call 'Room Off'. This turns the hearing aids' microphones off temporarily so you can only hear your mobile – allowing you to concentrate on the conversation without being disturbed by surrounding sounds. This can be easily deactivated again.

The simple and user-friendly RC-DEX remote control **gives you control of hearing aid functions** like volume and programs – at the simple click of a button. Extremely easy to use, the RC-DEX is small enough to fit on a key ring.

All our DEX devices are designed to be **intuitive and easy to use**. Particular attention has been paid to the user interface with simple, recognisable controls.

With DEX devices and CLEAR440-PA, a wireless world of near-natural sound is just a click away.



WIDEX CLEAR™ 440

STAY IN SYNC WITH LIFE

IT IS OUR MISSION, THROUGH
ORIGINALITY, PERSEVERANCE AND
RELIABILITY, TO DEVELOP **HIGH QUALITY**
HEARING AIDS THAT GIVE PEOPLE WITH A
HEARING LOSS THE SAME
OPPORTUNITIES FOR COMMUNICATION AS
THOSE WITH NORMAL HEARING.

WIDEX, WIDEX CLEAR440 and WidexLink
are trademarks of WIDEX A/S

www.widex.com/CLEAR440

WIDEX®

HIGH DEFINITION HEARING



Printed by KLS/2010-07
9 502 2135 001 #01



Mixed Sources
Product group from well-managed
forests and other controlled sources
www.fsc.org CertNo. SW-COC-003062
© 1996 Forest Stewardship Council