


PHONE-DEX

USER-FRIENDLY AND STYLISH



- 1. Talk button
- 2. End call and on/off button
- 3. Redial
- 4. Phone book

WWW.WIDEX.COM



Printed by FB/2011-07
9 502 2475 001 #01



PHONE-DEX

FOR EASY LANDLINE USE



A WHOLE NEW WORLD OF **WIRELESS** **SOUND**

The new DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.


They make it easy for your hearing aids to communicate with phones, TV's, remotes and other external devices.

PHONE-DEX

The PHONE-DEX is an all-purpose, cordless phone that streams crystal clear sound directly to your hearing aids. The PHONE-DEX provides you with good speech understanding in both ears.

It is incredibly easy to use - the only thing you need to do is hold the phone normally.

What's more, the PHONE-DEX works as a conventional phone as well. So family and friends can use it.



Getting started with PHONE-DEX couldn't be easier. It is ready to use as soon as you're ready to use it.

The PHONE-DEX can store up to twenty telephone numbers and comes with a base station that works as a charger. You can enjoy talking on the phone up to 300 metres away from the base station depending on surroundings.

