


T-DEX

FOR BLUETOOTH MOBILE PHONES



- 1. On/off, answering button
- 2. Volume button +/- (up/down)
- 3. Neck loop jack
- 4. Headphones jack
- 5. Bluetooth, battery and volume level indicator

WWW.WIDEX.COM





 Printed by FB/2012-01
 9 502 2498 001 #01



T-DEX

FOR BLUETOOTH MOBILE PHONES



A WHOLE NEW WORLD OF **WIRELESS** **SOUND**


The new DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.

They make communication between your hearing aids and your remote control, television, mobile or landline phone easy.

T-DEX

The T-DEX is a hands-free, wireless neck loop for use with all hearing aids with a telecoil. It lets you easily connect with Bluetooth mobile phones.

Simply place the T-DEX around your neck and then activate it when your phone rings. The sound is transmitted wirelessly to your hearing aids so you can hear and converse directly. Calls can be easily answered and ended by means of a single push button.



Because it amplifies the sound from your mobile phone, the T-DEX is perfect in situations where there is a lot of background noise. And its hands-free operation means that it is easy to talk on your phone while driving for example.

The T-DEX comes with a built-in rechargeable battery, as well as a charging cord.

