

EARMOULD/EAR-TIP FLEXIBILITY

Hearing loss profile	Hearing aid type	Receiver	Wire	Earmould/Ear-tip	Special features
Moderate to severe	S4-VS S2-VS	RIC P	EarWire RIC	Custom receiver soft ear-tip	Telecoil
				Custom receiver hard ear-tip	LED
				Instant receiver double ear-tip	Volume control
				Instant receiver tulip ear-tip	Program button
Severe to profound	S4-VS S2-VS	RITE SP	EarWire RITE	Custom hard earmould (Output Extender)	Size 675 battery
				Custom hard earmould (Regular sound bore)	DEX connectivity
		Wired RITE SP		Custom soft earmould (Regular sound bore)	Wireless fitting
					Classic program
					Audibility Extender
			WeatherCover		

WIDEX SUPER™

THE POWER TO HEAR

