

DEX™

A WIRELESS WORLD

T-DEX FOR BLUETOOTH MOBILE PHONES

The new DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.

The T-DEX is a hands-free, wireless neck loop for use with all hearing aids with a telecoil.

The sound is transmitted wirelessly to the hearing aids so the user can hear and converse directly. Calls can be easily answered and ended by means of a single push button.



Sound is transmitted wirelessly to MENU hearing aids from the T-DEX. It lets the user easily connect with Bluetooth mobile phones.

WIDEX[®]
HIGH DEFINITION HEARING

T-DEX
FOR BLUETOOTH MOBILE PHONES

T-DEX

FOR BLUETOOTH MOBILE PHONES



WIDEX[®]
HIGH DEFINITION HEARING

DEX[™]
A WIRELESS WORLD