

# WIDEX SUPER™

THE POWER TO HEAR



# WIDEX **SUPER**

## **Puts the sound where it should be**

### **- right in your ear**

WIDEX SUPER™ is a receiver-in-the-ear super power hearing aid for people with severe to profound hearing loss. We have designed it to be powerful but comfortable. It lets you hear speech clearly and naturally, and without distortion.

It can also help you determine where sounds are coming from\* and can even give you access to sounds that would otherwise be inaudible.

## **A real-world listening experience**

With WIDEX SUPER, you get a real-world listening experience. Being able to hear speech and conversation is of course essential, but in order to obtain a complete sound picture, you also need to hear everything from environmental sounds to background noise.

## **What you get**

A small, weather-resistant and stylish hearing aid that provides the advantages of comfort and clarity while giving you the power to hear.

# COLOURS



(074)  
Midnight black  
Charcoal grey



(071)  
Cappuccino brown  
Brown



(021)  
Warm beige  
Brown



(081)  
Tan silk  
Brown



(073)  
Titan grey  
Charcoal grey



(072)  
Winter silver  
Charcoal grey



(572)  
Winter silver  
Black

# HOW IT **SOUNDS**

For many people with severe to profound hearing loss, simply determining the direction of sound, if someone is approaching from the front or behind for example, is not possible.

Now, for the first time in a hearing aid of this type, SUPER gives you the possibility of identifying where sounds are coming from. It amplifies sounds in the range where your hearing is best. So the next time you're driving for example, you'll be able to better hear what people in the back of the car are saying.

## **Hear more**

More power does not always provide more hearing. That's why with SUPER, we have concentrated on audibility – giving you access to sounds that you may not have heard before. For example, SUPER features what we call the Audibility Extender. This is technology unique to Widex that actually moves high-pitched sounds such as birdsong and children's voices down to where you can hear them.





“My daughter’s voice sounds more natural and not so nasal as before.”

# WHAT IT **LOOKS** LIKE

New patent-  
pending  
WeatherCover  
(with nano-  
coating)

Volume  
control

Program button  
with LED





Telecoil

Receiver-  
in-the-ear  
(RITE)

Size 675  
battery

**WHERE** IT CAN HELP







### **Out and about**

If you are dependent on your hearing aids, then it's not good if you have to take them on and off all the time. SUPER is designed to be weather resistant, meaning you can use it all day long, no matter where or whatever the weather.

The microphones inside the hearing aid are also protected from moisture and wind noise.

### **On the phone**

The special Phone+ program allows you to easily hear phone calls directly through your hearing aids. Because your conversations are transmitted wirelessly, the sound quality is vastly improved.

### **In the classroom**

There is also a SUPER model with the option of an FM receiver so you can easily use it in places like classrooms or when listening to lectures.

### **All the time**

Even though it is small, SUPER has room for a long-life battery so you can enjoy more power for up to three weeks - without worrying about your hearing being interrupted.

# SUPER IS PERFECT IN MANY DIFFERENT LISTENING SITUATIONS

Listening situations	SUPER440	SUPER220
One-to-one	●●●●●	●●●●●
TV	●●●●●	●●●●●
Phone	●●●●●	●●●●○
Car	●●●●●	●●●●○
Outdoor	●●●●●	●●●○○
Classroom/lectures	●●●●●	●●●●●
Theatre/concerts	●●●●●	●●○○○
Meetings	●●●●●	●●●○○
Restaurant	●●●●●	●●○○○
Social events	●●●●●	●●○○○



VSD



VS

All SUPER hearing aids are compatible with our range of DEX assistive listening devices.



By choosing Widex hearing aids, you are choosing products from a company that has been WindMade certified. WindMade is the first global consumer label identifying companies that use wind power



Widex and SUPER are trademarks of Widex A/S.



Printed by FB/2012-09  
9 502 2929 001 #01

Widex A/S, Nymoellevej 6, DK-3540 Lyngø

[WWW.WIDEX.COM](http://WWW.WIDEX.COM)

**WIDEX**<sup>®</sup>  
HIGH DEFINITION HEARING