


WIDEX ZEN2GO™ TINNITUS DEVICE



The ultra-small WIDEX ZEN2GO tinnitus device is targeted at anyone suffering from tinnitus.

SUGGESTED FITTING RANGE



THE WIDEX ZEN2GO PACKAGE CONSISTS OF

- Two WIDEX ZEN2GO tinnitus devices, pre-fitted by the manufacturer to a flat 10 dB hearing loss
- Provided with S-receiver, earwire (size 2) and open ear-tip
- Matched RC-DEX
- Various cleaning accessories
- User instructions for WIDEX ZEN2GO and RC-DEX

WIDEX ZEN2GO FEATURES

- IE Zen
- SmartToggle
 - Zen Aqua
 - Zen Aqua + Noise
 - Zen Noise

WIDEX ZEN2GO OPTIONS

- Closed, vented or open fittings
- Instant and custom receiver ear-tips
- Can be fitted to provide additional amplification for people suffering from both tinnitus and hearing loss

USER FEEDBACK AND CONTROL

- SmartTone
- Battery alarm
- InterEar volume control (via RC-DEX)
- InterEar program shift (via RC-DEX)


ACCESSORIES

- RC-DEX included

C2-PA


MAXIMUM OUTPUT - EAR SIMULATOR

IEC 60118-0




MAXIMUM OUTPUT - 2CC COUPLER

IEC 60118-7 / ANSI S3.22-2009




OUTPUT - EAR SIMULATOR

IEC 60118-0



OUTPUT - 2CC COUPLER

IEC 60118-7 / ANSI S3.22-2009



Technical data

Typical data obtained through standard pure tone measurements.
Hearing aid set to Compass Test mode 1, unless stated otherwise.
Measured using a standard ITE coupler.

		IEC 60118-0	ANSI S3.22-2009 / IEC 60118-7
OSPL90	2500 Hz	114 dB SPL	104 dB SPL
	Peak	118 dB SPL	107 dB SPL
	Average	110 dB SPL	103 dB SPL
Acoustic output (Input 60 dB SPL)	2500 Hz	93 dB SPL	83 dB SPL
	Peak	94 dB SPL	84 dB SPL
	Average	85 dB SPL	78 dB SPL
Max gain (Input 50 dB SPL, Compass Max gain test mode)	2500 Hz	61 dB	52 dB
	Peak	64 dB	53 dB
	Average	58 dB	50 dB
Frequency range		100 Hz - 10000 Hz	100 Hz - 7800 Hz
Harmonic distortion (Input 70 dB SPL)	500 Hz	0.5%	0.3%
	800 Hz	0.5%	0.3%
	1600 Hz	0.7%	0.5%
Equivalent input noise		20 dB SPL	23 dB SPL
Battery drain (stand by)		1.13 mA	1.13 mA
Battery drain		1.14 mA	1.14 mA
Battery life / hours (Type 10 Zn-Air, 90 mAh)		80	80
Mobile phone immunity		IRIL: -33/-38 dB SPL	U-rating: M3

