

WIDEX **DRY-GO UV™**

ELECTRONIC DRYING STATION



WIDEX DRY-GO UV **- for perfect drying and cleaning**

The WIDEX DRY-GO UV™ electronic drying station is the perfect way to care for your hearing aids and earmoulds.

Small, portable and easy to use, the DRY-GO UV lets you clean and dry your hearing aids overnight.

And because it uses ultraviolet (UV) light, it is completely hygienic. So the risk of ear infections is reduced.

- Easy to use: start simply by touching the sensor button
- Switches off automatically after three hours of drying
- Energy efficient: uses very little power
- Standard charger: means that you can use it anywhere
- Hygienic: cleans with UV light
- Safe: UV light is switched off if the lid is opened
- Protects: extends the lifetime of your hearing aid







Perfect drying by convective flow of warm air



5V Micro-USB standard provides mobility



Hygienic cleansing by UV-C light for your health



Easy handling by touch of sensor button



Automatic room temperature modulation



Shuts down automatically after 3h drying



Energy-efficient and safe

WWW.WIDEX.COM



Printed by KLS/2013-08
9 502 3571 001 #01



MIX
Paper from
responsible sources
FSC® C022933



WIDEX[®]
HIGH DEFINITION HEARING