


# UNI-DEX

SIMPLE AND EASY TO USE



1. Push button/answer button
2. Room Off button
3. Neck loop
4. Microphone
5. Charging indicator
6. Charger connector
7. Jack plug



By choosing Widex hearing aids, you are choosing products from a company that has been WindMade certified. WindMade is the first global consumer label identifying companies that use wind power.

[WWW.WIDEXUSA.COM](http://WWW.WIDEXUSA.COM)



Widex and DEX are trademarks of Widex A/S.



Printed by CC/CBR173/0314  
9 502 3613 041 #01

Widex A/S, Nymoellevej 6, DK-3540 Lyngø

**WIDEX**<sup>®</sup>  
HIGH DEFINITION HEARING

# UNI-DEX

FOR PLUG AND PLAY MOBILE PHONE USE  
AND DIRECT AUDIO OUTPUT DEVICES



**DEX**<sup>™</sup>  
A WIRELESS WORLD

# A WHOLE NEW WORLD OF **WIRELESS** **SOUND**

The DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.

They make it easy to communicate between your hearing aids and audio equipment, remote controls, television, and mobile or landline phones.

## **UNI-DEX**

The UNI-DEX is a neck-worn **plug and play solution** for easy mobile phone use and other audio devices.

The built-in microphone in UNI-DEX allows you to **talk hands-free** when using a mobile device or Skype via PC or tablet.

And if the background noise becomes too annoying, the unique **Room Off feature** in UNI-DEX allows you to temporarily turn off the hearing aid microphones and concentrate on your conversation.

You can also use the UNI-DEX with your personal audio device.

The UNI-DEX is **easy to use**. Simply plug the audio jack into your mobile phone or audio device, hang the UNI-DEX around your neck, and you are ready to start listening. Streaming starts automatically.

Its long battery life means that UNI-DEX can stream for up to **40 hours**. And it takes only an hour to fully recharge again.

