

# EASYWEAR EARPIECE OVERVIEW

## FUSION, PASSION & SUPER

HEARING LOSS PROFILE	RECEIVER	SOLUTION TYPE	CUSTOM/INSTANT	EAR-TIP/MOULD TYPE	FAMILY/MODELS	MODELS
MILD	S	WIRED RIC	INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP	D-PA D-FS C-PA	PASSION FUSION
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP OPEN EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
MODERATE	M	WIRED RIC	INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP	D-PA D-FS C-PA C-FS	
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP OPEN EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
SEVERE	P	WIRED RIC	INSTANT	TULIP EAR-TIP DOUBLE EAR-TIP	D-FS C-FS VS/VSD	FUSION SUPER
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
	HP	RITE	CUSTOM	HARD EARMOULD**	D-FS C-FS	FUSION
		WIRED RITE	CUSTOM	SOFT EARMOULD		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
PROFOUND	SP	RITE	CUSTOM	HARD EARMOULD**	VS/VSD	SUPER
		WIRED RITE	CUSTOM	SOFT EARMOULD		
		SPECIAL*	CUSTOM	HARD EAR-TIP		

\*SPECIAL includes the modular and embedded types. In COMPASS GPS "SPECIAL" is used as acoustic identity for modular and embedded. \*\*Two options: Regular sound bore and Output Extender.

## INSTANT EAR-TIPS

S, M

S, M, P

### NEW

OPEN - SIZE: XS, S, L  
NEW SIZE XS

TULIP - SIZE: M, L  
NEW SIZE M

DOUBLE - SIZE: XS,S,M,L  
NEW SIZE XS



## CUSTOM EAR-TIPS/EARMOULD

S, M

S, M, P

S, M, P, HP, SP

OPEN HARD EAR-TIP

SOFT EAR-TIP

HARD EAR-TIP

HARD EAR-TIP/SPECIAL



## WIRED RECEIVERS

WIRED S, M AND P-RECEIVERS  
SIZE: 0, 1, 2, 3, 4  
LEFT/RIGHT



## RITE WIRES

## CUSTOM EARMOULD

SIZE: 0, 1, 2, 3, 4

RITE HP, SP

WIRED RITE HP, SP

HARD EARMOULD

SOFT EARMOULD



# EARPIECE OVERVIEW

## FUSION, PASSION, SUPER

**NEW!**

WIDEX **EASYWEAR™**

An extremely elegant, discreet receiver solution for FUSION, PASSION and SUPER