

WIDEX **CROS** FUSION



STANDARD TECHNOLOGY

- True-Input technology
- WidexLink - wireless connectivity
- Extremely long battery life
- Multiple EarWare options


The WIDEX CROS FUSION uses our unique wireless WidexLink technology and the new TRUE ISP, True-Input technology platform. Program button (transmission on/off) is standard. Uses a size 312 battery.

COMPATIBILITY TABLE



| | PASSION | FUSION | FASHION | 9 | <i>m</i> -CB | CIC | XP |
|----------|---------|--------|---------|---|--------------|-----|----|
| DREAM440 | x | x | x | x | x | x | x |
| DREAM330 | x | x | x | x | x | x | x |
| DREAM220 | x | x | x | x | x | x | x |
| DREAM110 | x | x | x | x | x | x | x |

WIDEX CROS COLOURS



KEY FEATURES

| | |
|--------------------------------------|------|
| Program button (transmission on/off) | X |
| Nano coating | X |
| Battery size | 312 |
| Water resistant | IP57 |
| Colours | 14 |

CROS-FS

PRODUCT DESCRIPTION:

| | |
|---|--------------|
| Frequency range | Up to 10 kHz |
| Battery size | 312 |
| Working current | ~ 1 mA |
| Battery life / hours (Type 312 Zn-Air, 145 mAh) | 150 (>135) |
| Microphone mode | Locator/Omni |
| Program button (transmission on/off) | ✓ |
| Nano coating | ✓* |
| Water resistant | IP57 |
| Input range | 113 dB SPL |

* Microphone cover

RADIO SPECIFICATIONS

| | |
|-------------------------|--------------------------|
| Frequency (range) | 10.6 MHz (10.2-11.0 MHz) |
| RF bandwidth | 660 kHz (-15 dB) |
| Channels | Single channel |
| Modulation | FSK |
| Radiated output power | 29pW (~ -75dBm) |
| Magnetic field strength | -54 dbuV/m @ 10 m |
| Duty cycle | 100 % Simplex |

| | |
|---|---|
| Conforms to the following standards/approvals | <ul style="list-style-type: none"> • EN/IEC 60950-1 • EN 300 330-2 • EN 301 489-1 • EN 301 489-3 • EN 62479 • FCC Part 15 C • FCC ID: TTY-DFS • CE Mark |
|---|---|

OVERVIEW

