

HEARING AID PERFORMANCE & FEATURE OVERVIEW

GLOBAL EDITION

AUDIBILITY/LISTENING IN QUIET		WIRELESS														NON-WIRELESS										
		UNIQUE440	UNIQUE330	UNIQUE220	UNIQUE110	DREAM440	DREAM330	DREAM220	DREAM110	DREAM100	DREAM50	DAILY100	DAILY50	CLEAR440	CLEAR330	CLEAR220	SUPER440	SUPER220	MENU10	MENU5	MENU3	MIND440	MIND330	MIND220	VITAL	BABY440
Variable Speed Compression	Automatic gearing for smooth sound and stable audibility	Yes (InterEar)	Yes (InterEar)	Yes (InterEar)	Yes																					
InterEar TruSound Compression	Featuring InterEar TruSound compression to support the localisation of sounds in complex environments					Yes							Yes				Yes									
TruSound Stabilisers	Adaptive compression in the TruSound system for optimal sound quality and to ensure speech is always heard in changing listening situations					Yes	Yes	Yes	Basic	Basic	Basic	Basic	Yes	Yes	Yes	Yes	Yes	Yes	Basic	Basic	Basic	Yes	Yes	Basic	Basic	Yes
Digital Pinna	Helps locate and focus naturally on the voice in front or behind by mimicking the normal pinna effect	Yes	Yes			Yes	Yes						Yes	Yes		Yes										
Enhanced Dynamic Range Compression	For audibility of soft sounds, speech understanding and comfort in loud environments	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Available via fine tuning	Available via fine tuning	Available via fine tuning	Yes	Yes	Yes	Yes	Yes
Fitting rationale options	WFR: Widex fitting rationale; NAL: NAL-NL2; DSL: DSL v5.0a	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL	WFR/NAL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL	WFR/DSL****	WFR/DSL****	WFR/DSL****	WFR/DSL	WFR/DSL	WFR/DSL	WFR	WFR/DSL	
LISTENING IN NOISE																										
True-Input Technology	Impressive speech intelligibility and sound quality in every situation due to an expanded input range of 113 dB SPL	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes													
High Definition Locator (Automatic Adaptive Directionality)	For better speech intelligibility in noise with channel-specific directionality for speech emphasis	15 ch	10 ch	6 ch	Broadband	15 ch	10 ch	5 ch	Broadband	Broadband	Broadband	Broadband	Broadband	15 ch	10 ch	5 ch	15 ch	5 ch	Add-on	Add-on	Add-on	15 ch	10 ch	5 ch	Omnidirectional	15 ch
Speech Enhancer RT	Reduces noise and enhances speech in real time in response to individual sound situations and hearing loss	Yes (InterEar)																								
Speech Enhancer SII	For understanding speech in complex environments; InterEar (IE) ensures speech enhancement for the dominant voice in a crowd		Yes			InterEar	Yes						InterEar	Yes		InterEar		Add-on****	Add-on****	Add-on****	Yes					Yes
Noise Reduction SIS	Provides listening comfort by reducing noise and intensifying speech	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Minimal	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Basic	Yes
Soft Level Noise Reduction	Unwanted background noises are reduced for greater listening comfort	Yes	Yes	Yes	Yes																					
FreeFocus (via M-DEX)	Enables the user to selectively focus on speech from the left side, right side or from behind	via COM-DEX APP	via COM-DEX APP	via COM-DEX APP	via COM-DEX APP	Yes	Yes	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes						
Reverse focus program	Choose listening to speech from behind	Yes	Yes	Yes		Yes	Yes	Yes					Yes	Yes	Yes	Yes	Yes									
Microphone Mode	Hearing aids can focus front and back for specific situations	Yes	Yes	Yes																						
USER COMFORT																										
Preference Control	User can control hearing aids to meet specific needs; either more audibility or more comfort	Yes	Yes	Yes	Yes																					
Sound Class Technology	Hearing aids adapt automatically to different listening situations for clarity, comfort and audibility	9 (InterEar)	5 (InterEar)	3	1																					
Multi-directional Active Feedback Cancelling	Minimises whistling; enhanced identification and cancelling of feedback	InterEar	InterEar	InterEar	Yes	InterEar	InterEar	InterEar	Yes	Yes	Yes	Yes	InterEar	InterEar	InterEar	InterEar	InterEar	InterEar	Yes	Yes	Yes	Yes	Yes	Yes	Static	Yes
TruSound Softener	Preserves comfort for sudden sounds at any input level by taking the edge off impact/impulse sounds	Yes	Yes	Yes		Yes	Yes						Yes	Yes		Yes										
TruSound AOC	Ensures distortion-free sound performance even when the hearing aid is fitted close to the fitting limit or high-level level output sounds	Yes	Yes	Yes	Yes	Yes	Yes	Yes					Yes	Yes	Yes	Yes	Yes	Yes				Yes	Yes			Yes
Hammer Detector	Detects and minimises undesirable loud transient sounds	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Wind Noise Attenuation	Speech and other sounds audible in wind	Yes																								
Wind noise control	Reduction of wind noise	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
USER CONTROLS																										
Listening programs	A wider selection of programs for more listening situations	5	4	3	3	5	4	3	2	3	3	3	3	5	4	3	5	3	2 + 3 Add-on	2 + 3 Add-on	2 + 3 Add-on	5	4	3	2	3
InterEar volume/program shift	Change the volume/program in both hearing aids simultaneously	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes								
Phone+	Hear phone calls in both hearing aids when they use WidexLink for transmission of sound	Yes	Yes			Yes	Yes						Yes	Yes		Yes	Yes									
SmartSpeak/SmartTone	Voice or tone alerts for easy hearing aid management	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Smart Tone only	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Add-on*	Add-on*	Add-on*	Yes	Yes	Yes	SmartTone	
Personal Zen	Flexible noise style adjustable according to desired frequency	Yes	Yes	Yes	Yes	Yes																				
Zen/Zen+	A tinnitus management tool for reducing the effects of tinnitus using soothing tones and/or noise; additional programs with Zen+	InterEar	InterEar	InterEar	Yes	InterEar	InterEar	InterEar	Yes	ZEN		ZEN	ZEN	InterEar/Yes	InterEar/Yes	Add-on	InterEar/Yes	Add-on	Add-on; Zen only	Add-on; Zen only	Add-on; Zen only	Yes	Add-on	Add-on		
Wireless CROS compatible	Wireless device to assist with single-sided deafness; utilises WidexLink					Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes													
HEARING HIGH FREQUENCIES																										
High-Frequency Boost	Clear and crisp sound in high frequencies (in music for example)	Yes																								
Personal Audibility Extender	High frequency sounds are audible in all listening situations and programs, using linear frequency transposition	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes																	
Audibility Extender	High frequency sounds are audible in one listening program, using linear frequency transposition	Yes	Yes	Yes	Yes								Yes	Yes	Yes	Yes	Yes	Yes	Add-on	Add-on	Add-on	Yes	Yes	Yes		Yes
FITTING AND FINE TUNING																										
Processing channels	For accurate processing and a more precise sound experience and speech intelligibility	15	10	6	4	15	10	5	3	3	3	3	3	15	10	5	15	5	10	5	3	15	10	5	2	15
Fitting frequency channels	Number of adjustable frequency bands	15	10	6	4	15	10	4	4	4	4	4	4	15	10	5	15	5	4	4	4	15	10	5	4	15
High-resolution Sensogram	Start the fitting process more accurately using in-situ thresholds for precision fitting	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Basic	Basic	Basic	Basic	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Sound Diary/Data logging	Allows easy monitoring of user data and hearing aid use - including characteristics of the listening environments and program selection	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Add-on	Add-on	Add-on	Yes	Yes	Yes		Yes
DEX DEVICES, CROS SOLUTIONS AND REMOTE CONTROLS																										
COM-DEX	A hands-free device for streaming high-quality sound from any mobile phone to hearing aids. COM-DEX APP available for iPhone and Android via COM-DEX	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes		Yes	Yes	Yes	Yes	Yes								
CALL-DEX	An ultra-compact device for streaming conversations from mobile phones to hearing aids	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes														
RC-DEX	Compact remote control for changing programs and volume	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes**	Yes	Yes	Yes	Yes	Yes	Yes	Yes								
UNI-DEX	Neck-worn plug & play solution for easy mobile phone use. Streams audio from any device with a mini-jack	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes								
TV-DEX	Improved TV listening without echo or delay	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Yes	Yes								
PHONE-DEX***	Seamless hearing with this landline telephone	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Yes	Yes								
FM+DEX	3 functions: FM, telecoil and audio line in	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes								
M-DEX	For better hearing on the mobile phone or other Bluetooth-compatible devices					Yes	Yes	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Yes	Yes								
T-DEX	Easy and affordable solution for connecting hearing aids with mobile phones. Connects via telecoil.	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes		Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
RC4-1	Up to 3 programs available for selection. Simultaneous volume adjustment for binaural fittings																		Yes	Yes	Yes	Yes	Yes	Yes		Yes
RC4-2	Up to 5 programs available for selection. Allows simultaneous or separate volume adjustment for binaural fittings																		Yes	Yes	Yes	Yes	Yes	Yes		Yes

PRODUCT OVERVIEW



* = SmartTone is included. SmartSpeak is an add-on. **Quickmatch only. *** PHONE-DEX only available in Europe and US/Canada ****ME-SP: only available with WFR ***** Not available for ME-SP



