

EARPIECE OVERVIEW

FUSION, PASSION AND SUPER

HEARING LOSS PROFILE	RECEIVER	SOLUTION TYPE	CUSTOM/INSTANT	EAR-TIP/MOULD TYPE	FAMILY/MODELS	MODELS
MILD	S	WIRED RIC	INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP	U-PA U-FS D-PA D-FS	PASSION FUSION
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP OPEN EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
MODERATE	M	WIRED RIC	INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP	U-PA U-FS D-PA D-FS	PASSION FUSION
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP OPEN EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
SEVERE	P	WIRED RIC	INSTANT	TULIP EAR-TIP DOUBLE EAR-TIP	U-FS D-FS S-VS/VSD	FUSION SUPER
			CUSTOM	SOFT EAR-TIP HARD EAR-TIP		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
	HP	RITE	CUSTOM	HARD EARMOULD**	U-FS D-FS	FUSION
		WIRED RITE	CUSTOM	SOFT EARMOULD		
		SPECIAL*	CUSTOM	HARD EAR-TIP		
PROFOUND	SP	RITE	CUSTOM	HARD EARMOULD**	S-VS/VSD	SUPER
		WIRED RITE	CUSTOM	SOFT EARMOULD		
		SPECIAL*	CUSTOM	HARD EAR-TIP		

FASHION, FASHION MINI, *m* AND 9-MODEL


HEARING LOSS PROFILE	TUBE	SOLUTION TYPE	CUSTOM/INSTANT	EAR-TIP/MOULD TYPE	FAMILY	MODELS
MILD TO SEVERE	THIN TUBE 0.9 ^{new}	THIN TUBE	CUSTOM	OPEN HARD EAR-TIP FLEX HARD EAR-TIP FLEX SOFT EAR-TIP	U-FA U-FM D-FA D- <i>m</i> CB	FASHION FASHION MINI <i>m</i>
	THIN TUBE 1.4 ^{new}		INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP		
	ÉLAN TUBE	THIN TUBE	CUSTOM	OPEN HARD EAR-TIP FLEX HARD EAR-TIP FLEX SOFT EAR-TIP	D-9	9
			INSTANT	OPEN EAR-TIP TULIP EAR-TIP DOUBLE EAR-TIP		
	CLASSIC TUBE		CUSTOM	OPEN HARD CLASSIC MOULD FLEX HARD CLASSIC MOULD FLEX SOFT CLASSIC MOULD	U-FA U-FM D-FA D- <i>m</i> CB D-9	FASHION FASHION MINI <i>m</i> 9

*SPECIAL includes the modular and embedded types. In compass GPS SPECIAL is used as acoustic identity for modular and embedded. **Two options : Regular sound bore and Output Extender.

EASYWEAR INSTANT EAR-TIPS



CUSTOM EAR-TIPS/EARMOULDS



EASYWEAR WIRED RECEIVERS



EASYWEAR THIN TUBES



RITE WIRES

CUSTOM EARMOULDS



EARPIECE OVERVIEW

FUSION, PASSION, SUPER,
FASHION, FASHION MINI, 9 AND *m*

NEW!
EASYWEAR THIN TUBE 0.9
AND
EASYWEAR THIN TUBE 1.4