

COM-DEX

KEEPS YOU CONNECTED. ALWAYS



A true hands-free device for streaming high-quality sound from any mobile phone to hearing aids. Can be combined with the COM-DEX app for convenient control of hearing aid volume, sound direction and program settings.



FEATURES

- True hands-free communication
- Compatible with any mobile phone
- Stream audio from all Bluetooth devices
- Available in three stylish colours
- 8 hours streaming time
- 8 days standby/control time*
- Comfortable fabric neck loop
- Control your hearing aids via COM-DEX app

* 12 hours per day

COM-DEX REMOTE MIC

HEAR CONVERSATIONS IN DETAIL



COM-DEX Remote Mic is a stylish add-on to your COM-DEX solution and is the answer to difficult one-to-one situations where clarity of speech is essential.

FEATURES

- Connects wirelessly via COM-DEX to all Widex wireless hearing aids
- Convenient steel clip for easy attaching
- Mute/Unmute microphone directly on the COM-DEX Remote Mic button
- Microphone can also be controlled using the COM-DEX app
- Up to 8 hours streaming time

TV-DEX

FOR ENJOYING TV



A palm-sized, user-friendly wireless assistive listening device purpose-built for enjoying TV and audio. Allows up to 10 hours of uninterrupted stereo TV or hi-fi sound between recharges.

FEATURES

- High quality audio
- Ultra-low delay "EchoFree™"
- Inputs for both TV and hifi systems
- Room Off feature
- 10 hours non-stop capacity
- Rechargeable using the TV-Base

PHONE-DEX

FOR EASY LANDLINE USE*



An all-purpose, cordless phone that streams crystal clear sound directly to your hearing aids. Provides the best possible listening experience for the user – without using a telecoil program or streamer. Also works as a conventional phone.

FEATURES

- Range to base up to 300 metres
- Phonebook for 20 numbers
- 100 hours of standby time and 10 hours call time on battery
- Can also be used as a regular phone – a phone for the whole family

*Available in EU and North America only

CALL-DEX

FOR DISCREET STREAMING FROM MOBILE PHONES



An ultra-compact device for streaming conversations from mobile phones to hearing aids. CALL-DEX simply inserts into jack plug of the phone and is then ready to use. Can stream up to 80 hours and is compatible with most mobile phones.

FEATURES

- Connects to all mobile phones with 3.5 mm jack output
- Streams automatically
- No buttons
- Constant streaming up to 80 hours
- Standby up to 3 weeks
- Small and easy to handle

UNI-DEX

FOR PLUG AND PLAY MOBILE PHONE USE



The UNI-DEX is a neck-worn plug and play solution for easy mobile phone use. The in-built microphone in UNI-DEX allows you to talk hands-free when using a mobile or Skype via PC or tablet.

FEATURES

- Connects to all devices with 3.5 mm jack output
- Streams automatically
- In-built microphone for hands-free talking
- Battery life 40 hours when streaming
- Battery recharges in 1 hour
- Room Off feature

FM+DEX

FOR FLEXIBLE STREAMING



A high quality streaming device that is small and discreet, and designed specifically for wireless hearing aids. Has three functions: FM, telecoil and audio line in, giving the user great flexibility. Easily operated with a toggle button.

FEATURES

- A flexible device for FM streaming
- Small and easy to handle
- Antenna range up to 30 m
- 10 hours battery capacity

RC-DEX

YOUR REMOTE CONTROL



A stylish, compact and user-friendly remote control that gives users better control of selected hearing aid functions. Allows the user to adjust volume and change between programs quickly and discreetly. Small enough to fit on a key ring. A lock switch prevents the remote being activated accidentally.

FEATURES

- Program toggle
- Volume up/down
- Key ring attachment
- Approximately 12 months of battery life
- Lock switch

DEX OVERVIEW

COMPATIBILITY	COM-DEX**	CALL-DEX	UNI-DEX	FM+DEX	RC-DEX	TV-DEX	PHONE-DEX
UNIQUE*	●	●	●	●	●	●	●
DREAM*	●	●	●	●	●	●	●
CLEAR	●		●	●	●	●	●
SUPER	●		●	●	●	●	●
WidexLink	●	●	●	●	●	●	●
Telecoil				●			
Bluetooth	●						
TV						●	
Audio	●		●	●		●	
Mobile phone	●	●	●				
Landline phone							●
Tamper resistant				●			
Compatible with other FM systems				●			

* not compatible with CIC micro. ** COM-DEX Remote Mic works only with COM-DEX



DEX™
A WIRELESS WORLD

OVERVIEW

A WHOLE WORLD OF WIRELESS SOUND



THE DEX FAMILY

The DEX™ communication accessories from Widex keep hearing aid users connected to a **whole world of sounds**. They make it easy for them to connect quickly and instantly to remote controls, TVs, personal audio systems, tablets and more.

Whether it is music, phone calls, TV or remote controls, there is **a DEX for every situation**.



WIDEX®
HIGH DEFINITION HEARING

WWW.WIDEX.COM

Widex A/S, Nymoellevvej 6, DK-3540 Lyngby

Google Play is a trademark of Google Inc.

mark of Apple Inc.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service

Widex and DEX are trademarks of Widex A/S.


Printed by KLS/2016-07. 9 502 4165 001 #01



9 502 416 500101



MIX
Paper from
responsible sources
FSC® C022833
www.fsc.org



KEEP CONNECTED WITH WIDEXLINK™

At the heart of the DEX communication accessories is our unique wireless technology, WidexLink.

It allows instant communication between hearing aids, direct connection to external audio devices, and coordination and synchronisation of data between two hearing aids.

WidexLink is designed specifically for the demands of hearing aids and DEX devices. It uses less power, is extremely fast and ensures there are no delays or dropout of sound.

