


POWER BTE Behind-the-ear

15/10E/10/5E/5/3E/3



Fitting range



Moderate to severe-to-profound hearing losses

Key features	15	10E	10	5E	5	3E	3
Processing channels	15	10	10	5	5	3	3
Linked Feedback Cancellation	Linked	Linked	Linked	Linked	•	•	•
Automatic Adaptive Directionality	15 ch.	10 ch.	10 ch.	5 ch.	Broadband	Broadband	Broadband
Real-Input	•	•	•	•	•	•	•
Noise Reduction w. Speech Focus	•	•	•	•	•	Basic	Basic
Listening programs	5	4	3	3	3	3	3
Linked Volume/Program Control	•	•	•	•	•	•	•
EasyVoice/EasyTone	•	•	•	•	•	EasyTone only	EasyTone only
Zen/Zen+	Linked	Linked	Linked	Linked	ZEN (no Zen+)		
Sound Stabiliser	•	•	•	•	Basic	Basic	Basic
Wireless CROS compatible	•	•	•	•	•	•	
RearView	•	•	•	•			
Linear frequency transposition	•	•	•	•			
Real Ear effect	•	•					
Impulse Noise Rejection II	•	•					
Linked Speech Enhancement System	Linked	•					
Phone Stream	•	•					
Linked Compression	•						
Personal Zen	•						
Accessories							
App for iPhone, Android and new functionality	Via COM-DEX	Via COM-DEX	Via COM-DEX	Via COM-DEX	Via COM-DEX		
DEX accessories	•	•	•	•	•	RC-DEX only	

Colours

- 021
 Warm Beige
- 081
 Tan Silk
- 072
 Winter Silver
- 073
 Titan Grey
- 074
 Midnight Black
- 071
 Capp. Brown
- 069
 Summer Gold

POWER BTE Behind-the-ear


15/10E/10/5E/5/3E/3

Technical data


		IEC 118-0	ANSI S3.22	
Output sound pressure level (Input 90 dB SPL)	2500 Hz	130	122	dB SPL
	Peak	138	134	dB SPL
	Average	133	125	dB SPL
Acoustic output (Input 60 dB SPL)	2500 Hz	116	108	dB SPL
	Peak	122	114	dB SPL
	Average	114	109	dB SPL
Full-on gain (Input 50 dB SPL)	2500 Hz	65	58	dB
	Peak	75	73	dB
	Average	70	62	dB
Telecoil output (Input 10 mA/m)	2500 Hz	106	99	dB SPL
	Peak	112	106	dB SPL
	Average	105	100	dB SPL
Frequency range		100 - 6250	100 - 6100	Hz
Total harmonic distortion (Input 70 dB SPL)	500 Hz	<2%	<2%	%
	800 Hz	<2%	<2%	%
	1600 Hz	<2%	<2%	%
Equivalent input noise level		20	20	dB SPL
Battery current (stand by)		0,93	0,93	mA
Battery current		1,01	1,03	mA
Battery life (Type 312 Zn-Air, 145 mAh)		265 (>225)	260 (>220)	h
Mobile phone immunity		IRIL: -5/12/9 dB SPL	U-rating: M4*/T4*	

*Low band


MAXIMUM OUTPUT - EAR SIMULATOR
IEC 118-0



MAXIMUM OUTPUT - 2CC COUPLER
ANSI S3.22 (2009)



OUTPUT - EAR SIMULATOR
IEC 118-0



OUTPUT - 2CC COUPLER
ANSI S3.22 (2009)

