

COM-DEX

KEEPS YOU CONNECTED. ALWAYS



A true hands-free device for streaming high-quality sound from any mobile phone to hearing aids. Can be combined with the COM-DEX app for convenient control of hearing aid volume, sound direction and program settings.



FEATURES

- True hands-free communication
- Compatible with any mobile phone
- Stream audio from all Bluetooth devices
- Available in three stylish colours
- 8 hours of streaming time
- 8 days of standby/control time*
- Comfortable fabric neck-loop
- Control your hearing aids via COM-DEX app

* 12 hours per day

COM-DEX REMOTE MIC

HEAR CONVERSATIONS IN DETAIL



COM-DEX Remote Mic is a stylish add-on to your COM-DEX solution and is the answer to difficult one-to-one situations, where clarity of speech is essential.

FEATURES

- Connects wirelessly via COM-DEX to all Widex wireless hearing aids
- Convenient steel clip for easy attaching
- Mute/Unmute microphone directly on the COM-DEX Remote Mic button
- Microphone can also be controlled using the COM-DEX app
- Up to 8 hours of streaming time

WIDEX TV PLAY™

FOR TV SOUND LIKE NEVER BEFORE



Outstanding sound performance in a discreet design. Easy to set up and compatible with WIDEX 2.4 GHz hearing aids. Provides direct, stable stereo streaming into the ears. Just play and stay connected.

FEATURES

- Never miss out on important TV moments
- Maximum streaming stability
- Outstanding sound
- Award-winning design

TV-DEX

FOR ENJOYING TV



A palm-sized, user-friendly wireless assistive listening device purpose-built for enjoying TV and audio. Allows up to 10 hours of uninterrupted stereo TV or hi-fi sound between recharges.

FEATURES

- High-quality audio
- Ultra-low delay "EchoFree™"
- Inputs for both TV and hi-fi systems
- Room Off feature
- 10 hours of non-stop capacity
- Rechargeable using the TV-Base

CALL-DEX

FOR DISCREET STREAMING FROM MOBILE PHONES



An ultra-compact device for streaming conversations from mobile phones to hearing aids. CALL-DEX simply inserted into the phone's jack plug and is then ready to use. Can stream up to 80 hours and is compatible with most mobile phones.

FEATURES

- Connects to all mobile phones with 3.5 mm jack output
- Streams automatically
- No buttons
- Constant streaming up to 80 hours
- Standby up to 3 weeks
- Small and easy to handle

UNI-DEX

FOR PLUG AND PLAY MOBILE PHONE USE



The UNI-DEX is a neck-worn plug and play solution for easy mobile phone use. The inbuilt microphone in UNI-DEX allows users to talk hands-free when using a mobile or Skype via PC or tablet.

FEATURES

- Connects to all devices with 3.5 mm jack output
- Streams automatically
- Inbuilt microphone for hands-free talking
- Battery life 40 hours when streaming
- Battery recharges in 1 hour
- Room Off feature

FM+DEX

FOR FLEXIBLE STREAMING



A high-quality streaming device that is small and discreet, and designed specifically for wireless hearing aids. Has three functions: FM, telecoil and audio line in, giving the users great flexibility. Easily operated with a toggle button.

FEATURES

- A flexible device for FM streaming
- Small and easy to handle
- Antenna range up to 30 m
- 10 hours of battery capacity

RC-DEX

YOUR REMOTE CONTROL



A stylish, compact and user-friendly remote control that gives users better control of selected hearing aid functions. Allows the users to adjust volume and change between programs quickly and discreetly. Small enough to fit on a keyring. A lock switch prevents the remote control being activated accidentally.

FEATURES

- Program toggle
- Volume up/down
- Keyring attachment
- Approximately 12 months of battery life
- Lock switch

ACCESSORIES OVERVIEW

COMPATIBILITY	TV PLAY***	TV-DEX	COM-DEX/ COM-DEX REMOTE MIC**	CALL-DEX	UNI-DEX	FM+DEX	RC-DEX
EVOKE*	●	●	●	●	●	●	●
BEYOND/BEYOND Z		●	●	●	●	●	●
UNIQUE*		●	●	●	●	●	●
DREAM*		●	●	●	●	●	●
SUPER		●	●		●	●	●
WidexLink		●	●	●	●	●	●
Telecoil						●	
Bluetooth			●				
TV	●	●					
Audio		●	●		●	●	
Mobile phone			●	●	●		
Landline phone							
Tamper-resistant						●	
Compatible with other FM systems						●	

* not compatible with CIC micro. ** COM-DEX Remote Mic works only with COM-DEX. *** EVOKE F2 only



A WHOLE WORLD OF WIRELESS SOUND OVERVIEW ACCESSORIES

WIDEX ACCESSORIES

Communication accessories from Widex keep hearing aid users connected to a **whole world of sound**.

They make it easy for users to connect quickly and instantly to remote controls, TVs, personal audio systems, tablets and more.

Be it music, phone calls, TV or remote controls, there is a **Widex accessory for every situation**.



WWW.WIDEX.COM


Widex A/S, Nymoellevvej 6, DK-3540 Lyngø

Google Play is a trademark of Google Inc.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Widex, DEX and TV PLAY are trademarks of WIDEX A/S.

Printed by ROS/2018-09, 9 502 4891 001 #01



KEEP CONNECTED WITH WIDEXLINK™

At the heart of our communication accessories is our unique wireless technology, WidexLink.

This technology allows instant communication between hearing aids, direct connection to external audio devices, and coordination and synchronisation of data between two hearing aids.

WidexLink is designed specifically for the demands of hearing aids and works with most Widex accessories. It uses less power, is extremely fast and ensures there are no delays or dropout of sound.

