



THIS SOUND
CHANGES
EVERYTHING

WIDEX MOMENT™

THIS SOUND CHANGES EVERYTHING





HEARING LIKE YOU REMEMBER

Even a great hearing aid has always sounded like a hearing aid – until now. The revolutionary WIDEX **MOMENT™** changes the game to deliver the most pure, natural sound ever.

Typically, sound processed in a hearing aid reaches the eardrum later than sound heard directly. And when these two ‘out of sync’ signals mix, you get an artificial sound. But now with WIDEX **MOMENT™** you can hear every moment like you used to.

PEOPLE
LOVE THE
SOUND



Are you tired of hearing aid 'revolutions' that no one can actually hear? We are. That's why WIDEX **MOMENT™** makes a real difference to your hearing in real life.

Tests show that a large majority of both new and experienced users prefer the Widex **PureSound™** program in WIDEX **MOMENT™** across different sound scenarios. Add this to the ability to effortlessly personalise your sound, and you have a hearing aid that will give you pure and natural sound like you've never heard before.





STREAM FROM **ALL YOUR DEVICES**

We've upgraded the 2.4 GHz Bluetooth® connectivity in WIDEX **MOMENT™** so that you can now stream from iOS, WIDEX **TV PLAY™*** and other **DEX** accessories, plus future connections with Android devices. So, now you're prepared for all the streaming entertainment the world has to offer in the highest quality possible.

* Available with WIDEX MOMENT™ compatibility late 2020.

BUILT TO LAST

WIDEX **MOMENT™** is designed for maximum moisture resistance and long life – since less corrosion means more durability.

You'll welcome advanced water-resistant nano coating to resist moisture. When there's no battery door, there are also fewer gaps for moisture to sneak into.





RECHARGEABILITY
IS NOW **DISCREET**

Convenience is discreet with WIDEX **MOMENT™**. Although the mRIC R D model is the market's smallest rechargeable receiver-in-canal hearing aid, you can rest easy knowing you'll have all the power you need to hear every moment.

No struggling with hearing aid batteries. Just a slim, simple charger that fits your life.



FINGERTIP CONTROL OF YOUR HEARING

With WIDEX **MOMENT™**, you have the option to fully control your hearing aid using the WIDEX **MOMENT™** app on your Apple or Android smartphone*. It's never been easier to fine-tune your hearing, so you get exactly the sound you prefer - in exactly the moments that are important to you.



WIDEX MOMENT™ App



Rigorous testing has verified that the app is easy and intuitive. And a 'dark mode' ensures discreet use.

* The WIDEX MOMENT™ app works with the following WIDEX MOMENT™ models: mRIC R D (MRR2D) and RIC 312 D (MRB2D). All other models use the WIDEX TONELINK™ app for control of key functions.

Menu
Reset sound
Frequency volume
In-app help



Change or select program
Mute volume and volume controls

WIDEX TONELINK™ App

GET MORE FROM YOUR WIDEX **MOMENT™**

WIRELESS ACCESSORIES – CONNECT, CONTROL AND COMMUNICATE

DEX accessories help you effortlessly connect to audio equipment, remote controls, television and mobile or landline phones. Ask your hearing care professional for more details.



COM-DEX
For hands-free
streaming



**COM-DEX
REMOTE MIC**
For improved
audibility in noise



RC-DEX
Your remote
control



TV-DEX
For enjoying TV



FM-DEX
For flexible
streaming



red dot award 2018
winner



**STYLISH AND STABLE TV STREAMING
WITH WIDEX TV PLAY™***

* Available with WIDEX MOMENT™ compatibility late 2020.

WIDEX **MOMENT™** MODELS AND STYLES

Widex **MOMENT™** is available in three behind-the-ear (RIC) models and three in-the-ear (ITE) models. Your hearing care professional can help you choose the ideal model for your specific hearing loss and lifestyle.

BEHIND-THE-EAR (RIC) MODELS



RIC 10



RIC 312 D



mRIC R D

IN-THE-EAR (ITE) MODELS



XP



CUSTOM ITE (Im/IP)



CIC

CHOOSE YOUR **COLOURS**


























Choose between 13 different colours for RIC models and three standard colours for ITE models.



NEW COLOURS!
BASED ON A GLOBAL USER
HAIR COLOUR STUDY

WHY WIDEX **MOMENT™** IS RIGHT FOR YOU

To help you hear your best, WIDEX **MOMENT™** features exclusive Sound Class Technology that automatically adapts your hearing aid to your specific listening requirements. Without pushing any buttons, WIDEX **MOMENT™** ensures optimal clarity, comfort and audibility in each different situation.

| PERFORMANCE LEVEL/BENEFIT | SOCIAL | | MUSIC | | QUIET | | PARTY | | TRANSPORT | | URBAN | |
|---------------------------|---|---|---|---|---|--|---|---|---|---|--|-------------|
| | | | | Classical or contemporary | | with speech | | with speech | | with speech | | with speech |
| WIDEX MOMENT™ 440 |  |  |  |  |  |  |  |  |  |  |  | |
| WIDEX MOMENT™ 330 | |  |  | |  |  |  |  | |  | | |
| WIDEX MOMENT™ 220 | |  |  | | |  | |  | | | | |
| WIDEX MOMENT™ 110 | |  |  | | |  | | | | | | |

| PERFORMANCE LEVEL/BENEFIT | Improved speech understanding Get accentuated speech and reduced noise from the surroundings | Wind noise reduction Hear clearly even when it's windy | Natural sound experience Choose between the renowned Widex sound or the new Widex PureSound™ | Automatic sound adjustment No need to adjust the hearing aid in every situation | Direct Bluetooth® streaming Stream sound directly from Apple or Android devices* | Easy to control Intuitive app and Artificial Intelligence guide you to the best sound | WIDEX REMOTE CARE Enables fitting and fine-tuning at any location the client requests | Real-life insights Facilitates counselling by providing access to the clients' personal programs |
|----------------------------------|--|--|--|---|--|---|---|--|
| WIDEX MOMENT™ 440 | | | | | YES | YES | YES | YES |
| WIDEX MOMENT™ 330 | | | | | YES | YES | YES | YES |
| WIDEX MOMENT™ 220 | | | | | YES | YES | YES | YES |
| WIDEX MOMENT™ 110 | | | | | YES ** | YES ** | YES | YES |

* Prepared for future connections with Android (ASHA)

** Doesn't include **SoundSense Learn**



Printed 2020-01. 9 502 5335 001 #01

Widex, WIDEX MOMENT™, Puresound™, ZeroDelay™ Technology, TruAcoustics™, WIDEX REMOTE CARE™, ToneLink™, TV Play™ and DEX are trademarks of Widex A/S.



WIDEX A/S, Nymøllevej 6, DK-3540 Lyngø

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android is a trademark of Google Inc. Google Play and the Google Play logo are trademarks of Google Inc.