



By choosing Widex hearing aids, you are choosing products from a company that has been WindMade certified. WindMade is the first global consumer label identifying companies that use wind power.



Widex and DEX are trademarks of Widex A/S.



Printed by FB/2012-03
9 502 2415 001 #01

Widex A/S, Nymoellevej 6, DK-3540 Lyngø

WWW.WIDEX.COM

WIDEX[®]
HIGH DEFINITION HEARING

FM+DEX

COMPACT AND EASY TO USE



1. FM
2. Audio
3. Telecoil
4. Volume
5. On/Off
6. Toggle



FM+DEX

FM+DEX

FOR FLEXIBLE STREAMING



DEXTM

A WIRELESS WORLD



A WHOLE NEW WORLD OF **WIRELESS** **SOUND**

The new DEX™ assistive listening devices from Widex let you take advantage of a whole new world of sounds.

They make communication between your hearing aids and your remote control, television, mobile or landline phone, or FM system easy.

FM+DEX

The FM+DEX is a high-quality streaming device designed specifically for hearing aids.

It converts signals from an FM receiver, telecoil or line in and streams them to your Widex wireless hearing aids.

It is compact and user-friendly, and you can easily operate it with a toggle button for program selection and volume control.

A rechargeable battery in the FM+DEX provides you with up to ten hours of streaming.

It is also designed with an optional function for locking program and volume and comes with a retention solution suitable for children.

The FM+DEX is ideal in a variety of situations such as the classroom, meetings, lectures and concerts.



