

RC-DEX


QUICK, DISCREET AND HELPFUL



* Available also in black

- 1. LED
- 2. Volume up
- 3. Program toggle
- 4. Volume down

WWW.WIDEXUSA.COM



Printed by FB/2010-12
9 502 1838 041 #02

CBR124/1210



TAKE CONTROL OF YOUR HEARING AIDS!

RC-DEX

The easy-to-use, keychain-sized wireless remote that puts you in charge



A NEW WORLD OF **WIRELESS** **SOUND**

3 simple push-button controls on a lightweight device that is small enough to fit on a keychain!
That's the RC-DEX!

With the RC-DEX, you can have *control* over your hearing aids – and all completely *wirelessly*. The unique technology at the heart of our RC-DEX means that your hearing aids can communicate easily and instantly with your remote control.

The RC-DEX is so smart, it can even send signals to your hearing aids and speak to you to inform you of volume and program change.

RC-DEX

The RC-DEX is a user-friendly and very small hearing aid remote control. Its simple and intuitive design gives you better control of hearing aid features.

An LED indicator on the front of the remote control shows when it is ON. The RC-DEX even gives you discreet control if you use it from *inside* your pocket or purse!

